

# Imaging Beyond Consumer Cameras – Proseminar (911.422)

## Exercise sheet C

### Image registration

#### Exercise 1.

20 P.

In this multi-step exercise, you will learn how to use image registration to *automatically segment an MR image of the human brain into medically relevant regions*. In medical image analysis, the pipeline that is outlined below is referred to as *atlas-based image segmentation*.

**Background.** The idea is the following: say, we would have an MR image, where medical experts have hand-segmented certain regions (e.g., the brain stem, the ventricles, etc.). In medical imaging, such an image is often referred to as an *atlas* (see figure below), or a *template*. Sometimes, the atlases are constructed from images of *multiple* patients, or from a *single* patient only.



To automatically segment a new image (from a new patient), e.g., to study the size of certain brain structures related to some disease (e.g., Alzheimer's, etc.), we could try (1) to *map* this patients' MR image onto the atlas, then (2) transfer the labels from the atlas to the new image, and eventually transform the MR image back into its original space.

**Caution:** An important detail that we need to take care of is to *strip* the skull from the MR images, since image registration would immediately hinge onto the skull structure, instead of the soft-tissue structures (because of very strong gradients). For affine registration this is not a big issue, however, when a more advanced registration model (e.g., non-linear) is used, *skull-stripping* typically makes a huge difference.

#### The exercise is split into multiple parts:

1. Download example brain MRI scans (CUMC12 dataset) from here [🔗](#) and select 1-2 images
2. Download the SPL 2008 brain atlas from here [🔗](#)
3. Register the image **affinely** to the atlas (using `aladdin` of NiftyReg for instance)
4. Register the image **non-linearly** to the atlas (using `f3d` of NiftyReg for instance)
5. Visualize the registration result (just one slice) via 3DSlicer's *Checkerboard Filter*
6. Transform labels from atlas to image (using `reg_resample` of NiftyReg for instance)
7. Extract the *corpus callosum* (using `Convert3D` (`c3d`); read the documentation for how to extract a label with a certain number)
8. Visualize the corpus callosum (e.g., using ParaView, or ITKSnap)
9. Document the process in 1-2 pages with visualizations and parameter settings of the tools that you used

## Resources

Below, you can find links to software packages that you can use (I do recommend using NiftyReg in combination with BET):

- [!\[\]\(95b42f0077faf7439a26242a54e021ec\_img.jpg\) NiftyReg](#)
- [!\[\]\(e097ab4c08b8186dd0908330bbc2dc28\_img.jpg\) ANTs](#)
- [!\[\]\(1e9d865c5de095f8e3304757c49e79d7\_img.jpg\) 3DSlicer](#)
- [!\[\]\(735b10d724a5f0ec5005c4eb3eb9c9d1\_img.jpg\) FSL Brain Extraction Tool \(BET\)](#)

## Additional help

Note that the registration and label transfer steps can be (almost) exactly replicated from the [Segmentation Propagation Tutorial](#).

Some of the images come in the format of a .img and a .hdr file. In case NiftyReg can't directly process these, I do recommend converting them using Convert3D (c3d), e.g., using

```
c3d OAS1_0001_MR1_mpr_n4_anon_111_t88_masked_gfc.img OAS1_0001_MR1_mpr_n4_anon_111_t88_masked_gfc.nii
```

I would also check if the conversion has worked by (1) visually checking the registration results and (2) checking the meta information, e.g., via the `-info` parameter of Convert3D.

**Evaluation criteria.** 5 points for affine registration; +5 points for non-linear registration; +5 points for transforming the labels back and +5 points for visualization.